



145 James Way  
 Southampton, Pennsylvania 18966 USA  
**P** 215-526-2151 • 888-755-4370  
**F** 215-526-2167 • 888-258-4293

**ESOP** Employee Owned  
 for Your Benefit

[www.advantapure.com](http://www.advantapure.com)

PURITY IN FLUID FLOW SYSTEMS®

## Change Notification

Dear Valued Customer:

This letter serves as formal notification of changes that may affect products or services purchased by your organization.

### 1. General Information

Date Issued:	March 26, 2019
Subject:	Nordson Medical Changes to the PharmLok Sanitary Clamps (TC Style) and PharmaLok Tube Clamps. Nordson Change Notification # PCN390
Change Notification Number:	1932
Effective Dates:	March through July 2019

### 2. Change Description

Part Series Affected:	CPAX150-U71, CPAX150-U75, CPAX150-U706 CPAN75-U71, CPAN75-U75, CPAN75-U706 CPT100-U71 and CPT90-U71
Description of Changes:	Minor changes to PharmaLok Sanitary Clamps (TC Style) CPAX150 & CPAN 75 are being made to the hinge area to improve the operation of the clamp when opened to its furthest position. Changes to PharmaLok Tube Clamps CPT100 & CT90 are being made in the latch area to prevent the upper half of the clamp from swiping to the side of the lower clamp when pinching down onto any tubing.
Old Specifications Were:	See Nordson Medical Change Notification
New Specifications Are:	See Nordson Medical Change Notification

### 3. Summary Checklist: Impact to Customer

<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No	Drawing Change?
<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No	Material Change?
<input checked="" type="checkbox"/> Yes	<input type="checkbox"/> No	Form Change? Minimal Impact on Form or Fit
<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No	Documentation Change?
<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No	Certification Change?
<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No	Packaging/Labeling Change?
<input checked="" type="checkbox"/> Yes	<input type="checkbox"/> No	Other? Process Improvement – Minor Changes - will improve the Clamps function.
<input type="checkbox"/>	<input type="checkbox"/>	

Should you have any questions regarding this notice, please contact your Customer Service Representative at 215-526- 2300. Other change notifications and details about our Change Notification Policy can be found on our website at <http://www.advantapure.com/change-notifications.htm>. Thank you for your continued support of NewAge® Industries AdvantaPure.

Alex Kakad, Product Marketing Manager

Signature and Title

QC-FORM-0002, Revision E, 2/26/2018

March 26, 2019

Date